


## Anti-Fadd Monoclonal Antibody For Wb, Ihc-P, Elisa - Q61160

Item Number: CM0002075



### Introduction

Explore the role of FADD in apoptosis and innate immunity with this validated rabbit monoclonal antibody, ideal for WB, IHC-P, and ELISA across human, mouse, and rat samples.

[Learn More](#)

| Parameter        | Value                                                                                                           |
|------------------|-----------------------------------------------------------------------------------------------------------------|
| Product Name     | FADD Rabbit mAb                                                                                                 |
| Remarks/Alias    | Mort1/FADD; FADD                                                                                                |
| Species          | Mouse                                                                                                           |
| GeneID (Human)   | 8772                                                                                                            |
| GeneID           | 14082                                                                                                           |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 1-205 of mouse FADD (NP_034305.1) |
| Source           | Rabbit                                                                                                          |
| Category         | Monoclonal Antibodies                                                                                           |
| Application      | WB, IHC-P, ELISA                                                                                                |
| Cross Reactivity | Human, Mouse, Rat                                                                                               |
| SWISS            | Q61160                                                                                                          |
| Protein Weight   | 23 kDa                                                                                                          |
| Shipping         | Ice bag                                                                                                         |