


Anti-Fadd Rabbit Polyclonal Antibody For Wb, Ip, Elisa - Q13158

Item Number: CM0028429



Introduction

Rabbit polyclonal antibody against human FADD (FAS-associated death domain protein).

Validated for WB, IP, and ELISA; cross-reacts with mouse and rat FADD. Useful for apoptosis, death receptor signaling, and innate immunity research.

[Learn More](#)

Specification	Details
Product Name	FADD Rabbit pAb
Remarks/Alias	GIG3; IMD90; MORT1; FADD
Species	Human
Gene ID (Human)	8772
Gene ID	8772
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-208 of human FADD (NP_003815.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13158
Protein Weight	23kDa
Shipping	Ice bag