

Anti-Fads1 Rabbit Polyclonal Antibody For Wb, Elisa - O60427

Item Number: CM0023754

Introduction

Rabbit polyclonal antibody against human FADS1 (delta-5 desaturase) for WB and ELISA. Validated in human, mouse, rat. Involved in PUFA metabolism and inflammatory response; ideal for lipid metabolism research.

[Learn More](#)

Parameter	Value
Product Name	FADS1 Rabbit pAb
Remarks/Alias	D5D, TU12, FADS6, FADSD5, LLCDL1
Species	Human
GeneID (Human)	3992
GeneID	3992
Immunogen	Synthetic peptide within amino acids 100-200 of human FADS1 (NP_037534.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O60427
Protein Weight	52kDa
Shipping	Ice bag