


Anti-Fak (Focal Adhesion Kinase) Rabbit Monoclonal Antibody For Wb, If/Icc, Elisa - Q05397

Item Number: CM0003534



Introduction

High-quality rabbit monoclonal antibody targeting human FAK (Focal Adhesion Kinase 1), validated for WB, IF/ICC, ELISA. Reacts with human, mouse, rat. Ideal for cell adhesion and migration research.

[Learn More](#)

Parameter	Value
Product Name	FAK Rabbit mAb
Remarks/Alias	FAK; FADK; FAK1; FRNK; FADK 1; PPP1R71; p125FAK; pp125FAK
Species	Human
GeneID (Human)	5747
GeneID	5747
Immunogen	Synthetic peptide
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q05397
Protein Weight	119kDa
Shipping	Ice bag