


Anti-Fak Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Ip, Elisa - Q05397

Item Number: CM0016968



Introduction

Rabbit polyclonal antibody against human FAK (Q05397) validated for WB, IHC-P, IF/ICC, IP, ELISA. Cross-reacts with human, mouse, rat. Suitable for cell adhesion, migration, and cancer research.

[Learn More](#)

Parameter	Value
Product Name	[KD Validated] FAK Rabbit pAb
Remarks/Alias	FAK; FADK; FAK1; FRNK; FADK 1; PPP1R71; p125FAK; pp125FAK
Species	Human
GeneID (Human)	5747
GeneID	5747
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 734-834 of human FAK (NP_722560.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q05397
Protein Weight	119kDa
Shipping	Ice bag