


Anti-Fbp2 Rabbit Polyclonal Antibody For Wb, If-P, Elisa - O00757

Item Number: CM0024795



Introduction

FBP2 rabbit polyclonal antibody validated for WB, IF-P, ELISA; reacts with human, mouse, rat. Suitable for studying fructose-1,6-bisphosphatase isozyme 2 in glycogen synthesis and muscle metabolism.

[Learn More](#)

Parameter	Value
Product Name	FBP2 Rabbit pAb
Remarks/Alias	CORLK; FBP2
Species	Human
Gene ID (Human)	8789
Gene ID	8789
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-339 of human FBP2 (NP_003828.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Mouse, Rat
SWISS	O00757
Protein Weight	37kDa
Shipping	Ice bag