

Anti-Fbx32/Fbxo32 Polyclonal Antibody For Wb, If-P, Elisa - Q969P5

Item Number: CM0025101

Introduction

Rabbit polyclonal antibody against human Fbx32/FBXO32 (Atrogin-1, MAFbx). Validated for WB, IF-P, ELISA. Cross-reacts with human and mouse. Targets an E3 ubiquitin ligase involved in muscle atrophy.

[Learn More](#)

Parameter	Value
Product Name	Fbx32/FBXO32 Rabbit pAb
Remarks/Alias	Fbx32; MAFbx; Fbx32/FBXO32
Species	Human
Gene ID (Human)	114907
Gene ID	114907
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 206-355 of human Fbx32/FBXO32 (NP_478136.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q969P5
Protein Weight	42kDa
Shipping	Ice bag