


Anti-Fbxo7 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q9Y3I1

Item Number: CM0027177



Introduction

Rabbit polyclonal antibody detects FBXO7 in human, mouse, and rat. Suitable for WB, IF/ICC, and ELISA. FBXO7 is an F-box protein involved in ubiquitination and mitophagy, with a predicted molecular weight of 59kDa.

[Learn More](#)

Parameter	Value
Product Name	FBXO7 Rabbit pAb
Remarks/Alias	FBX, FBX7, PKPS, FBXO7, PARK15, FBXO7
Species	Human
GeneID (Human)	25793
GeneID	25793
Immunogen	Recombinant fusion protein containing amino acids 353-522 of human FBXO7 (NP_036311.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Y3I1
Protein Weight	59kDa
Shipping	Ice bag