


Anti-FcεR1G Polyclonal Antibody For Wb, Elisa - P30273

Item Number: CM0013730



Introduction

Rabbit polyclonal antibody against FCER1G (Fc receptor gamma-chain) for Western blot and ELISA. Reacts with human, mouse, and rat samples. Suitable for allergy and immunology research.

[Learn More](#)

Parameter	Value
Product Name	FCER1G Rabbit pAb
Remarks/Alias	FCRG; FCER1G
Species	Human
GeneID (Human)	2207
GeneID	2207
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 19-86 of human FCER1G (NP_004097.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P30273
Protein Weight	10kDa
Shipping	Ice bag