

Anti-Fcgr2B Rabbit Polyclonal Antibody For Wb, If-P, Elisa - P31994

Item Number: CM0016422



Introduction

Rabbit polyclonal antibody targeting human FCGR2B (CD32), validated for WB, IF-P, and ELISA. Detects human, mouse, and rat. Essential for studying FcγRIIB-mediated immune regulation and autoimmune disease mechanisms.

[Learn More](#)

Parameter	Value
Product Name	FCGR2B Rabbit pAb
Remarks/Alias	CD32; FCG2; CD32B; FCGR2; IGFR2; FCGR2C; FcGRIIB; FcRII-c; FcγgammaRIIb; FCGR2B
Species	Human
GeneID (human)	2213
GeneID	2213
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 43-217 of human FCGR2B (NP_003992.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P31994
Protein Weight	34kDa
Shipping	Ice bag