

Anti-Fcgr3B Rabbit Polyclonal Antibody For Wb, Elisa - 075015

Item Number: CM0016212



Introduction

FCGR3B (CD16b) rabbit polyclonal antibody validated for WB and ELISA. Recognizes human and rat low-affinity IgG Fc receptor III-B, a neutrophil-expressed immune complex trap. Sourced from rabbit, immunized with recombinant human FCGR3B protein (1-233aa).

[Learn More](#)

Core Product Specifications	Details
Product Name	FCGR3B Rabbit pAb
Remarks/Alias	CD16; FCG3; CD16A; CD16b; FCGR3; CD16-I; FCGR3A; FCR-10; FCRIII; FCRIIIb; FCGR3B
Species	Human
GeneID (Human)	2215
GeneID	2215
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-233 of human FCGR3B (NP_000561.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Rat
SWISS	O75015
Protein Weight	26kDa
Shipping	Ice bag