

Anti-Fermt2 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q96Ac1

Item Number: CM0025654



Introduction

Rabbit polyclonal antibody against human FERMT2 (Kindlin-2, MIG-2). Validated for WB, IHC-P, ELISA. Reacts with human, mouse, rat. Suitable for studying integrin signaling, focal adhesion, and Wnt pathway.

[Learn More](#)

Parameter	Value
Product Name	FERMT2 Rabbit pAb
Remarks/Alias	MIG2; KIND2; mig-2; UNC112; PLEKHC1; UNC112B; FERMT2
Species	Human
GeneID (Human)	10979
GeneID	10979
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 380-550 of human FERMT2 (NP_006823.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q96AC1
Protein Weight	78kDa
Shipping	Ice bag