


## Anti-Ferritin Heavy Chain Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P02794

Item Number: CM0016187



### Introduction

CamelBio supplies this anti-Ferritin Heavy Chain Rabbit Polyclonal Antibody for WB, IHC-P, IF/ICC, and ELISA. Cross-reacts with human, mouse, rat FTH1 (21 kDa). Ideal for iron metabolism and neurodegeneration research.

[Learn More](#)

| Parameter        | Value                                                                                                                           |
|------------------|---------------------------------------------------------------------------------------------------------------------------------|
| Product Name     | Ferritin Heavy Chain Rabbit pAb                                                                                                 |
| Remarks/Alias    | FHC; FTH; HFE5; PLIF; FTHL6; PIG15; Ferritin Heavy Chain                                                                        |
| Species          | Human                                                                                                                           |
| GeneID (Human)   | 2495                                                                                                                            |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 1-183 of human Ferritin Heavy Chain (NP_002023.2) |
| Source           | Rabbit                                                                                                                          |
| Category         | Polyclonal Antibodies                                                                                                           |
| Application      | WB, IHC-P, IF/ICC, ELISA                                                                                                        |
| Cross Reactivity | Human, Mouse, Rat                                                                                                               |
| SWISS            | P02794                                                                                                                          |
| Protein Weight   | 21kDa                                                                                                                           |
| Shipping         | Ice bag                                                                                                                         |