


Anti-Fgf1 Monoclonal Antibody For Wb, If/Icc, Elisa - P05230

Item Number: CM0008800



Introduction

Rabbit monoclonal antibody targeting human FGF1 (aFGF, heparin-binding growth factor 1). Validated for WB, IF/ICC, ELISA applications. Cross-reacts with human, mouse, and rat. Ideal for investigating cell proliferation and angiogenesis pathways.

[Learn More](#)

Parameter	Value
Product Name	FGF1 Rabbit mAb
Remarks/Alias	AFGF; ECGF; FGFA; ECGFA; ECGFB; FGF-1; HBGF1; HBGF-1; GLIO703; ECGF-beta; FGF-alpha; FGF1
Species	Human
GeneID (Human)	2246
GeneID	2246
Immunogen	Synthetic Peptide (A synthetic peptide corresponding to a sequence within amino acids 56-155 of human FGF1 (P05230))
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P05230
Protein Weight	17kDa
Shipping	Ice bag