

Anti-Fgf1 Rabbit Polymab® For Wb, If/Icc, Ihc-P, Elisa - P05230

Item Number: CM0010912



Introduction

Recombinant rabbit monoclonal FGF1 antibody (PolymAb®) validated for WB, IF/ICC, IHC-P, and ELISA. Detects human, mouse, and rat FGF1 with a predicted MW of 17 kDa. Ideal for angiogenesis and cell signaling research.

[Learn More](#)

Parameter	Value
Product Name	FGF1 Rabbit PolymAb®
Remarks/Alias	AFGF; ECGF; FGFA; ECGFA; ECGFB; FGF-1; HBGF1; HBGF-1; GLIO703; ECGF-beta; FGF-alpha
Species	Human
GeneID (Human)	2246
GeneID	2246
Immunogen	Recombinant protein
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P05230
Protein Weight	17kDa
Shipping	Ice bag