

Anti-Fgf2 Rabbit Polyclonal Antibody For Wb, Elisa - P09038

Item Number: CM0022336



Introduction

FGF2 Rabbit Polyclonal Antibody validated for WB and ELISA. Detects human FGF2 and cross-reacts with mouse and rat. Suitable for research on cell growth, differentiation, and angiogenesis.

[Learn More](#)

Parameter	Value
Product Name	FGF2 Rabbit pAb
Remarks/Alias	BFGF; FGFB; FGF-2; HBGF-2; F2
Species	Human
GeneID (human)	2247
GeneID	2247
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 100-200 of human FGF2 (NP_001997.5).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P09038
Protein Weight	31kDa
Shipping	Ice bag