

Anti-Fggy Polyclonal Antibody For Wb, Elisa - Q96C11

Item Number: CM0013991

Introduction

Rabbit polyclonal antibody targeting human FGGY carbohydrate kinase domain-containing protein. Validated for WB and ELISA; cross-reacts with human, mouse, and rat. Suitable for D-ribulose phosphorylation and metabolic repair studies.

[Learn More](#)

Parameter	Value
Product Name	FGGY Rabbit pAb
Remarks/Alias	FGGY
Species	Human
GeneID (Human)	55277
GeneID	55277
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-250 of human FGGY (NP_060761.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q96C11
Protein Weight	60kDa
Shipping	Ice bag