


Anti-Fgr Polyclonal Antibody For Wb, Ihc-P, Elisa - P09769

Item Number: CM0019845



Introduction

Rabbit polyclonal antibody against human FGR (SWISS P09769) for Western blot, IHC-P, and ELISA. Cross-reacts with human and rat. Suitable for investigating non-receptor tyrosine kinase signaling in immune responses and cytoskeleton dynamics.

[Learn More](#)

Parameter	Value
Product Name	FGR Rabbit pAb
Remarks/Alias	SRC2; c-fgr; c-src2; p55-Fgr; p58-Fgr; p55c-fgr; p58c-fgr; FGR
Species	Human
GeneID (Human)	2268
GeneID	2268
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 1-100 of human FGR (NP_005239.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Rat
SWISS	P09769
Protein Weight	59kDa
Shipping	Ice bag