

Anti-Fibronectin Rabbit Monoclonal Antibody For Wb, If/Icc, Elisa - P02751

Item Number: CM0002015



Introduction

Fibronectin Rabbit Monoclonal Antibody (P02751) validated for WB, IF/ICC, ELISA. Reacts with human and mouse. Ideal for cell adhesion, wound healing, and extracellular matrix studies in diverse research applications.

[Learn More](#)

Parameter	Value
Product Name	Fibronectin Rabbit mAb
Remarks/Alias	FN; CIG; FNZ; MSF; ED-B; FINC; GFND; LETS; GFND2; SMDCF; Fibronectin
Species	Human
Gene ID	2335
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-2477 of human FN1 (P02751)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	P02751
Protein Weight	272kDa
Shipping	Ice bag