

Anti-Fkbp10 Rabbit Polymab® For Wb, If/Icc, Elisa - Q96Ay3

Item Number: CM0011652

Introduction

FKBP10 Rabbit PolymAb® is a recombinant rabbit monoclonal antibody that detects human, mouse, and rat FKBP10. Validated for WB, IF/ICC, and ELISA, it is ideal for studying protein folding, collagen biosynthesis, and bone development.

[Learn More](#)

Parameter	Value
Product Name	FKBP10 Rabbit PolymAb®
Remarks/Alias	OI6; BRKS; OI11; TLH1; BRKS1; FKBP65; PPIASE; hFKBP65
Species	Human
GeneID (Human)	60681
GeneID	60681
Immunogen	Recombinant protein
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q96AY3
Protein Weight	64kDa
Shipping	Ice bag