

Anti-Fkbp12 Rabbit Monoclonal Antibody For Wb And Elisa - P62942

Item Number: CM0003046



Introduction

High-quality rabbit monoclonal antibody against human FKBP12, validated for WB and ELISA. Specific for human, mouse, rat. Ideal for TGF-beta signaling, protein folding, and ryanodine receptor studies.

[Learn More](#)

Parameter	Value
Product Name	FKBP12 Rabbit mAb
Remarks/Alias	FKBP1; PKC12; PKC12; FKBP12; PPIASE; FKBP-12; FKBP-1A
Species	Human
Gene ID (Human)	2280
Gene ID	2280
Immunogen	Synthetic peptide within amino acids 29-108 of human FKBP12 (P62942)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P62942
Protein Weight	12kDa
Shipping	Ice bag