

Anti-Flad1 Rabbit Polyclonal Antibody For Wb, Elisa - Q8Nff5

Item Number: CM0013160



Introduction

Rabbit polyclonal antibody against human FLAD1 (FAD synthase/diphosphatase), validated for WB and ELISA. Detects endogenous levels of ~65 kDa FLAD1 protein in human and mouse samples.

[Learn More](#)

Parameter	Value
Product Name	FLAD1 Rabbit pAb
Remarks/Alias	FAD1; FADS; PP591; LSMFLAD; FLAD1
Species	Human
Gene ID	80308
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 388-587 of human FLAD1 (NP_079483.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q8NFF5
Protein Weight	65kDa
Shipping	Ice bag