

Anti-Flcn Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q8Nfg4

Item Number: CM0012938



Introduction

Rabbit polyclonal FLCN antibody validated for WB, IF/ICC, and ELISA; reacts with human, mouse, and rat. Ideal for studying mTORC1 signaling, autophagy, and Birt-Hogg-Dubé syndrome research.

[Learn More](#)

Parameter	Value
Product Name	FLCN Rabbit pAb
Remarks/Alias	BHD; FLCL; DENND8B; FLCN
Species	Human
GeneID	201163
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 341-579 of human FLCN (NP_659434.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8NFG4
Protein Weight	64kDa
Shipping	Ice bag