

Anti-Fndc4 Polyclonal Antibody For Wb, Elisa - Q9H6D8

Item Number: CM0018325

Introduction

This FNDC4 rabbit polyclonal antibody detects human, mouse, and rat FNDC4 and is validated for WB and ELISA. Useful for studying inflammation, insulin sensitivity, and adipose biology.

[Learn More](#)

Parameter	Value
Product Name	FNDC4 Rabbit pAb
Remarks/Alias	FRCP1; FNDC4
Species	Human
GeneID (Human)	64838
GeneID	64838
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 45-167 of human FNDC4 (NP_073734.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9H6D8
Protein Weight	25kDa
Shipping	Ice bag