

Anti-Fndc4 Rabbit Polyclonal Antibody For Ihc-P/Elisa - Q9H6D8

Item Number: CM0024588

Introduction

Rabbit polyclonal antibody against human FNDC4 (Fibronectin type III domain-containing protein 4), validated for IHC-P and ELISA. Cross-reacts with mouse. Ideal for inflammation, glucose tolerance, and adipose biology research.

[Learn More](#)

Parameter	Specification
Product Name	FNDC4 Rabbit pAb
Remarks/Alias	FRCP1; FNDC4
Species	Human
Gene ID	64838
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 45-167 of human FNDC4 (NP_073734.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9H6D8
Protein Weight	25kDa
Shipping	Ice bag