

Anti-Fntb Rabbit Polyclonal Antibody For Wb, Elisa - P49356

Item Number: CM0020667



Introduction

Rabbit polyclonal antibody against human FNTB, validated for WB and ELISA. Recognizes human, mouse, and rat FNTB, a farnesyltransferase critical for protein prenylation and signal transduction.

[Learn More](#)

Parameter	Details
Product Name	FNTB Rabbit pAb
Remarks/Alias	FPTB; FNTB
Species	Human
Gene ID	2342
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-120 of human FNTB (NP_002019.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P49356
Protein Weight	49kDa
Shipping	Ice bag