

Anti-Fpgs Monoclonal Antibody For Wb, Elisa - Q05932

Item Number: CM0003315

Introduction

Validated for WB and ELISA, this rabbit monoclonal antibody detects human and rat FPGS, the key enzyme in folate polyglutamation and methotrexate metabolism. Ideal for research on one-carbon metabolism, antifolate pharmacology, and mitochondrial enzymology.

[Learn More](#)

Parameter	Value
Product Name	FPGS Rabbit mAb
Species	Human
GeneID (Human)	2356
GeneID	2356
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 450-587 of human FPGS (Q05932)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Rat
SWISS	Q05932
Protein Weight	65kDa
Shipping	Ice bag