


Anti-Fshr Rabbit Polyclonal Antibody For Wb, If-P, Elisa - P23945

Item Number: CM0016141



Introduction

Rabbit polyclonal antibody targeting human FSHR (Follicle-stimulating hormone receptor). Validated for WB, IF-P, ELISA and cross-reacts with mouse and rat. Ideal for reproductive physiology and GPCR pathway studies.

[Learn More](#)

Parameter	Value
Product Name	FSHR Rabbit pAb
Remarks/Alias	LGR1; ODG1; FSHR1; FSHRO; FSHR
Species	Human
Gene ID	2492
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 631-695 of human FSHR (NP_000136.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P23945
Protein Weight	78kDa
Shipping	Ice bag