

Anti-Fsp1/S100A4 Rabbit Monoclonal Antibody For Wb, If/Icc, If-P, Ihc-P, Elisa - P26447

Item Number: CM0002504



Introduction

Rabbit monoclonal antibody targeting human S100A4 (FSP1), validated for Western blot, immunofluorescence, immunohistochemistry, and ELISA. Cross-reacts with mouse and rat. Ideal for studying metastasis, cell motility, and apoptosis.

[Learn More](#)

Parameter	Value
Product Name	FSP1/S100A4 Rabbit mAb
Remarks/Alias	42A; 18A2; CAPL; FSP1; MTS1; P9KA; PEL98; FSP1/S100A4
Species	Human
Gene ID (Human)	6275
Gene ID	6275
Immunogen	Recombinant protein corresponding to amino acids 2-101 of human S100A4 (P26447)
Source	Rabbit
Category	Monoclonal Antibody
Application	WB, IF/ICC, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P26447
Protein Weight	12kDa
Shipping	Ice bag