


Anti-Fubp1 Rabbit Monoclonal Antibody For Ip, If/Icc, If-P, Ihc-P, Elisa - Q96Ae4

Item Number: CM0005469



Introduction

Rabbit monoclonal antibody against human FUBP1, suitable for IP, IF/ICC, IF-P, IHC-P, and ELISA. Detects human, mouse, and rat FUBP1. Generated from a synthetic peptide corresponding to amino acids 545-644. For transcription regulation research.

[Learn More](#)

Parameter	Value
Product Name	FUBP1 Rabbit mAb
Remarks/Alias	FBP; FUBP; hDH V; FUBP1
Species	Human
GeneID (Human)	8880
GeneID	8880
Immunogen	Synthetic Peptide: A synthetic peptide corresponding to a sequence within amino acids 545-644 of human FUBP1 (Q96AE4).
Source	Rabbit
Category	Monoclonal Antibodies
Application	IP, IF/ICC, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q96AE4
Protein Weight	68kDa
Shipping	Ice bag