


Anti-Fut8 Polyclonal Antibody For Wb, Elisa - Q9Byc5

Item Number: CM0015023



Introduction

High-quality rabbit polyclonal antibody against human FUT8 (alpha-1,6-fucosyltransferase). Validated for WB and ELISA, cross-reacts with mouse. Ideal for glycobiology, IgG research, and CDG studies.

[Learn More](#)

Parameter	Value
Product Name	FUT8 Rabbit pAb
Remarks/Alias	CDGF, CDGF1, FUT8
Species	Human
Gene ID	2530 (Human)
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 376-575 of human FUT8 (NP_835367.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9BYC5
Protein Weight	67kDa
Shipping	Ice bag