


Anti-Fxn / Frataxin Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q16595

Item Number: CM0013948



Introduction

Rabbit polyclonal antibody against human frataxin (FXN). Validated for WB, IHC-P, IF/ICC, and ELISA. Reacts with human, mouse, and rat. Suitable for studies of iron-sulfur cluster biogenesis and Friedreich's ataxia.

[Learn More](#)

Parameter	Value
Product Name	FXN / Frataxin Rabbit pAb
Remarks/Alias	FA; X25; CyaY; FARR; FRDA; FXN / Frataxin
Species	Human
GeneID (Human)	2395
GeneID	2395
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 42-210 of human FXN / Frataxin (NP_000135.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q16595
Protein Weight	23kDa
Shipping	Ice bag