


Anti-Fxr1 Rabbit Polyclonal Antibody For Wb And Elisa - P51114

Item Number: CM0029943



Introduction

Rabbit polyclonal antibody against human FXR1, suitable for WB and ELISA with cross-reactivity to mouse and rat. Ideal for investigating mRNA translation regulation, neurogenesis, and muscle development.

[Learn More](#)

Parameter	Value
Product Name	FXR1 Rabbit pAb
Remarks/Alias	FXR1P; CMYP9A; CMYP9B; MYOPMIL; MYORIBF; FXR1
Species	Human
GeneID (Human)	8087
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 250-350 of human FXR1 (NP_005078.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P51114
Protein Weight	70kDa
Shipping	Ice bag