

Anti-Fxr2 Polyclonal Antibody For Wb, If/Icc, Elisa - P51116

Item Number: CM0025279



Introduction

Rabbit polyclonal antibody against human FXR2 (P51116), validated for WB, IF/ICC, and ELISA. Cross-reacts with mouse and rat. Suitable for mRNA stability and hippocampal neurogenesis research.

[Learn More](#)

Parameter	Value
Product Name	FXR2 Rabbit pAb
Remarks/Alias	FXR2P; FMR1L2; FXR2
Species	Human
Gene ID (Human)	9513
Gene ID	9513
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 550-650 of human FXR2 (NP_004851.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P51116
Protein Weight	74kDa
Shipping	Ice bag