

Anti-Gaba Transaminase (Abat) Rabbit Polyclonal Antibody For Wb, If-P, Elisa - P80404

Item Number: CM0030592

Introduction

Rabbit polyclonal antibody against human GABA transaminase (ABAT). Validated for WB, IF-P, and ELISA with cross-reactivity to mouse and rat. Ideal for studying GABA metabolism and neurological disorders.

[Learn More](#)

Parameter	Value
Product Name	GABA transaminase (ABAT) Rabbit pAb
Remarks/Alias	GABAT; NPD009; GABA-AT; GABA transaminase (ABAT)
Species	Human
GeneID (Human)	18
GeneID	18
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P80404
Protein Weight	56kDa
Shipping	Ice bag