


Anti-Gabarap Polyclonal Antibody For Wb, Ihc-P, Elisa - 095166

Item Number: CM0013774



Introduction

GABARAP Rabbit Polyclonal Antibody validated for WB, IHC-P, ELISA. Detects human, mouse, rat GABARAP, a key autophagy and GABA(A) receptor transport protein. Ideal for neurology and oncology research.

[Learn More](#)

Parameter	Value
Product Name	GABARAP Rabbit pAb
Remarks/Alias	MM46, ATG8A, GABARAP-a, GABARAP
Species	Human
Gene ID (Human)	11337
Gene ID	11337
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1-100 of human GABARAP (NP_009209.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O95166
Protein Weight	14kDa
Shipping	Ice bag