


Anti-Gadd45A Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P24522

Item Number: CM0021059



Introduction

Rabbit polyclonal antibody against GADD45A (DDIT1), validated for WB, IHC-P, and ELISA. Reacts with human, mouse, rat species. Suitable for DNA damage and cell cycle research applications.

[Learn More](#)

Parameter	Value
Product Name	GADD45A Rabbit pAb
Remarks/Alias	DDIT1; GADD45; GADD45A
Species	Human
GeneID (Human)	1647
GeneID	1647
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 1-165 of human GADD45A (NP_001915.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P24522
Protein Weight	18kDa
Shipping	Ice bag