

Anti-Gal3St4 Rabbit Polyclonal Antibody For Wb And Elisa - Q96Rp7

Item Number: CM0026059

Introduction

Rabbit polyclonal antibody targeting human GAL3ST4 (galactose-3-O-sulfotransferase 4) for Western blot and ELISA applications. Cross-reacts with human, mouse, and rat. Suitable for studying sulfation pathways in O-linked glycoproteins.

[Learn More](#)

Parameter	Value
Product Name	GAL3ST4 Rabbit pAb
Remarks/Alias	GAL3ST-4; GAL3ST4
Species	Human
Gene ID	79690
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 137-486 of human GAL3ST4 (NP_078913.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q96RP7
Protein Weight	54kDa
Shipping	Ice bag