

Anti-Galectin-1 (Lgals1) Monoclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P09382

Item Number: CM0009326



Introduction

Anti-Galectin-1 rabbit monoclonal antibody validated for WB, IHC-P, IF/ICC, ELISA in human, mouse, rat. Ideal for apoptosis, cell proliferation, and immunoregulation studies, including T-cell apoptosis and Th17 differentiation.

[Learn More](#)

Parameter	Value
Product Name	Galectin 1/LGALS1 Rabbit mAb
Remarks/Alias	GBP; GAL1; Galectin 1/LGALS1
Species	Human
GeneID (Human)	3956
GeneID	3956
Immunogen	A synthetic peptide within amino acids 1-100 of human Galectin-1 (P09382)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P09382
Protein Weight	15kDa
Shipping	Ice bag