

Anti-Galk1 Rabbit Polyclonal Antibody For Wb, Elisa - P51570

Item Number: CM0024492



Introduction

Rabbit polyclonal antibody against human GALK1, raised against recombinant protein (aa 1-120). Validated for WB and ELISA, cross-reacts with human and mouse. Ideal for galactose metabolism and kinase research.

[Learn More](#)

Parameter	Value
Product Name	GALK1 Rabbit pAb
Remarks/Alias	GK1; GALK; HEL-S-19; GALK1
Species	Human
GeneID (Human)	2584
GeneID	2584
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-120 of human GALK1 (NP_000145.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P51570
Protein Weight	42kDa
Shipping	Ice bag