


Anti-Gata4 Polyclonal Antibody For Wb, Ihc, If/Icc And Elisa - P43694

Item Number: CM0018012



Introduction

Rabbit polyclonal antibody against human GATA4, validated for WB, IF/ICC, IF-P, IHC-P, and ELISA. Reacts with human, mouse, and rat samples. Suitable for cardiac development and transcription factor studies.

[Learn More](#)

Parameter	Value
Product Name	GATA4 Rabbit pAb
Remarks/Alias	TOF; ASD2; VSD1; TACHD; GATA4
Species	Human
Gene ID	2626
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 341-440 of human GATA4 (NP_002043.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P43694
Protein Weight	45kDa
Shipping	Ice bag