


Anti-Gch1 Polyclonal Antibody For Wb, If/Icc, Elisa - P30793

Item Number: CM0017572



Introduction

Rabbit polyclonal antibody against human GCH1, validated for WB, IF/ICC, and ELISA. Detects native GCH1 in human, mouse, and rat samples. Ideal for tetrahydrobiopterin biosynthesis, dopamine synthesis, and pain sensitivity research.

[Learn More](#)

Parameter	Value
Product Name	GCH1 Rabbit pAb
Remarks/Alias	GCH; DYT5; DYT14; DYT5a; GTPCH1; HPABH4B; GTP-CH-1; GCH1
Species	Human
GeneID (Human)	2643
GeneID	2643
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 55-250 of human GCH1 (NP_000152.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P30793
Protein Weight	28kDa
Shipping	Ice bag