

Anti-Gclm Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P48507

Item Number: CM0027300



Introduction

Rabbit polyclonal antibody against human GCLM (glutamate-cysteine ligase modifier subunit), validated for WB, IF/ICC, and ELISA in human, mouse, and rat samples. Ideal for glutathione biosynthesis research.

[Learn More](#)

Parameter	Value
Product Name	GCLM Rabbit pAb
Remarks/Alias	GLCLR, GCLM
Species	Human
Gene ID (Human)	2730
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-274 of human GCLM (NP_002052.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P48507
Protein Weight	31kDa
Shipping	Ice bag