


Anti-Gcn2 Polyclonal Antibody (Rabbit) For Wb, Ihc-P, If/Icc, Elisa - Q9P2K8

Item Number: CM0016880



Introduction

High-quality rabbit polyclonal antibody against human GCN2 (EIF2AK4), validated for WB, IHC-P, IF/ICC, and ELISA. Cross-reacts with mouse and rat. Ideal for amino acid starvation and stress response research.

[Learn More](#)

Parameter	Value
Product Name	GCN2 Rabbit pAb
Remarks/Alias	GCN2; PVOD2
Species	Human
GeneID (Human)	440275
GeneID	440275
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9P2K8
Protein Weight	187kDa
Shipping	Ice bag