


Anti-Gcn2 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q9P2K8

Item Number: CM0026279



Introduction

Rabbit polyclonal antibody against human GCN2 (EIF2AK4), suitable for WB, IF/ICC, ELISA. Cross-reacts with mouse and rat. GCN2 is a stress-sensing kinase that phosphorylates eIF-2 α , regulating translation during amino acid starvation.

[Learn More](#)

Parameter	Value
Product Name	GCN2 Rabbit pAb
Remarks/Alias	GCN2; PVOD2
Species	Human
GeneID (Human)	440275
GeneID	440275
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-300 of human GCN2 (NP_001013725.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9P2K8
Protein Weight	187kDa
Shipping	Ice bag