


Anti-Gdf-8/Myostatin Polyclonal Antibody For Wb, Ihc-P, Elisa - O14793

Item Number: CM0029826



Introduction

High-quality rabbit polyclonal antibody against GDF-8/Myostatin, validated for Western blot, immunohistochemistry, and ELISA. Suitable for human, mouse, and rat samples. Ideal for muscle growth research.

[Learn More](#)

Parameter	Value
Product Name	GDF-8/Myostatin Rabbit pAb
Remarks/Alias	GDF8; MSLHP; GDF-8/Myostatin
Species	Human
Gene ID (Human)	2660
Gene ID	2660
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 267-375 of human GDF-8/Myostatin (NP_005250.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O14793
Protein Weight	43 kDa
Shipping	Ice bag