

Anti-Gdf15 Monoclonal Antibody For Wb And Elisa - Q99988

Item Number: CM0007137



Introduction

Rabbit monoclonal antibody targeting human growth/differentiation factor 15 (GDF15), validated for Western blot and ELISA. Ideal for research on stress response, metabolic regulation, cancer cachexia, and nausea mechanisms.

[Learn More](#)

Parameter	Value
Product Name	GDF15 Rabbit mAb
Remarks/Alias	PDF; MIC1; PLAB; MIC-1; NAG-1; PTGFB; GDF-15; GDF15
Species	Human
Gene ID (Human)	9518
Gene ID	9518
Immunogen	Synthetic peptide within amino acids 100-200 of human GDF15 (NP_004855.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q99988
Protein Weight	34kDa
Shipping	Ice bag