


Anti-Gif/Intrinsic Factor Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P27352

Item Number: CM0029743



Introduction

Rabbit polyclonal antibody against human GIF (Intrinsic Factor), validated for WB, IF/ICC, and ELISA. Cross-reacts with mouse. Suitable for cobalamin absorption and gastric mucosa research.

[Learn More](#)

Parameter	Value
Product Name	GIF Rabbit pAb
Remarks/Alias	IF; GIF; INF; IFMH; TCN3
Species	Human
Gene ID (Human)	2694
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 20-300 of human GIF (NP_005133.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	P27352
Protein Weight	45kDa
Shipping	Ice bag