

Anti-Gipr Polyclonal Antibody For Wb, Ihc-P, Elisa - P48546

Item Number: CM0027349



Introduction

High-quality rabbit polyclonal antibody against human GIPR (Gastric inhibitory polypeptide receptor), validated for WB, IHC-P, ELISA. Cross-reacts with mouse and rat. Ideal for GIP signaling research.

[Learn More](#)

Parameter	Value
Product Name	GIPR Rabbit pAb
Remarks/Alias	PGQTL2; GIPR
Species	Human
GeneID (human)	2696
GeneID	2696
Immunogen	Recombinant protein / Recombinant fusion protein containing a sequence corresponding to amino acids 22-138 of human GIPR (NP_000155.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P48546
Protein Weight	53kDa
Shipping	Ice bag