

## Anti-Gkn1 Polyclonal Antibody For Wb, Elisa - Q9Ns71

Item Number: CM0013403

### Introduction

Rabbit polyclonal antibody against human GKN1, validated for WB and ELISA with cross-reactivity to mouse. Gastroke-1 is a mitogenic protein involved in gastric mucosal integrity. Ideal for gastrointestinal research.

[Learn More](#)

| Parameter        | Value                                    |
|------------------|------------------------------------------|
| Product Name     | GKN1 Rabbit pAb                          |
| Remarks/Alias    | FOV; CA11; AMP18; BRICD1; foveolin; GKN1 |
| Species          | Human                                    |
| GeneID (Human)   | 56287                                    |
| GeneID           | 56287                                    |
| Immunogen        | Recombinant protein                      |
| Source           | Rabbit                                   |
| Category         | Polyclonal Antibodies                    |
| Application      | WB, ELISA                                |
| Cross Reactivity | Mouse                                    |
| SWISS            | Q9NS71                                   |
| Protein Weight   | 20kDa                                    |
| Shipping         | Ice bag                                  |